

TAKE A FRESH LOOK AT IMPLANTOLOGY

---



**Anthogyr**  
PRIME MOVER IN IMPLANTOLOGY

# Anthogyr

PRIME MOVER IN IMPLANTOLOGY



# axiom<sup>®</sup>

MULTI<sub>LEVEL</sub><sup>®</sup>



In perfect harmony and continuity with all that has been proposed until now, the **Axiom<sup>®</sup> Multi Level<sup>®</sup>** range opens up the extent of possibilities in offering **total compatibility** between the **Bone Level** and **Tissue Level** philosophies.

Audacious and innovative, yet simple and accessible, this new range offers a new freedom, an unprecedented approach to dental implantology with unparalleled and complete fluidity.



## NEW IMPLANT

Axiom® TL, Tissue Level

Platform diameters  
Ø4.0 - Ø4.8 mm,  
Gingival heights  
1.5 - 2.5 - 3.5 mm

inLink® connection for  
multiple screw-retained  
restorations

Tissue Favored Design

Indexed connection for  
single-unit restoration

REG or PX profile

### Biological Safety

Prosthetic management at the gingival level  
A unique epithelial attachment  
Tissue maximization: Tissue Favored Design anatomic profile

### Prosthetic Ease

Access and prosthetic manipulations made easier  
Time saved by carrying out a one-step surgery  
Large choice of gingival profiles

### Comfort of screw-retained prostheses

Enlarged indications thanks to Angulated Access channels  
No cement contamination  
Comfortable maintenance and scalability

axiom®  
MULTI LEVEL®

POWER IN YOUR HANDS

### SMART RANGE

Total integration into  
the Axiom® world

### FULL COMPATIBILITY

between Axiom® TL and Axiom® BL

Axiom® TL,  
Tissue Level

Axiom® BL,  
Bone Level

SERENITY of a Sameda®  
full CAD-CAM prosthesis

Sina Z® Zirconia, Titanium or  
Cobalt-Chromium prosthesis

inLink® abutment

## INNOVATION

New inLink® connection

Integrated  
inLink®  
lock system

Machined compartment  
for inLink® system

Fixation  
lock

Angulated  
Access

Retaining  
ring

### Ease

inLink® connection by locks integrated in the prosthesis

### Flexibility

Correction of implant axes divergences without intermediary abutment

### Comfort

Angulated Access Channels up to 25°

### Smart Range

Complete compatibility between Axiom® TL, Tissue Level and Axiom® BL, Bone Level

### Serenity

Adapted, fitting and reliable custom-made Sameda® prosthesis

Anthogyr  
PRIME MOVER IN IMPLANTOLOGY



# Anthogyr

PRIME MOVER IN IMPLANTOLOGY

## Axiom® Multi Level®

More than a new system, a new philosophy...  
More than a new implant, **your new power.**

Anthogyr  
2 237, avenue André Lasquin  
74700 Sallanches  
Phone +33 (0)4 50 58 02 37  
[www.anthogyr.com](http://www.anthogyr.com)



[www.axiom-multilevel.com](http://www.axiom-multilevel.com)